


Technical Specification

Product: SARS-CoV-2 Antigen Kit (Colloidal Gold) (for self-testing)



Product name	SARS-CoV-2 Antigen Kit (Colloidal Gold) (for self-testing)	Device type	In vitro diagnostic medical device
Intended use	For qualitative detection of SARS-CoV-2 nucleocapsid antigen to aid in the rapid diagnosis of SARS-CoV-2 infections	Test principle	Rapid immunochromatographic assay
Pack size	1 test/box, 5 tests/box, 25 tests/box	Catalog number	CG123001, CG123005, CG123025
Sample	Nasal swab	Time to result	12 minutes
Qualitative / Quantitative	Qualitative	Reader required	None. Visually read
Automation	No, manual	Storage condition	2-30°C, away from light
Kit components	Test cartridge Sterilized swab Extraction tube (prefilled with buffer) Instructions for Use	Shelf-life	18 months from manufacturing date

Product photo

1T

REF CG123001



5T

REF CG123005



Packaging Information

Cat No.	Picture	Description	Pack Size (Test)
CG123001		SARS-CoV-2 Antigen Kit (Colloidal Gold) 480 single kits\carton 66*44*34cm net weight:14KG	1 test per kit
CG123005		SARS-CoV-2 Antigen Kit (Colloidal Gold) 1000 single kits\carton 67*38*37cm net weight:15KG	5 tests per kit
CG123025		SARS-CoV-2 Antigen Kit (Colloidal Gold) 1000 tests\carton 77*36*39cm, net weight:15KG	25 tests per kit



CERTIFICATE

EC Certificate No. 1434-IVDD-483/2021

**EC Design-examination
Directive 98/79/EC concerning
in vitro diagnostic medical devices**

Polish Centre for Testing and Certification certifies
that manufactured by:

Goldsite Diagnostics Inc.

**No. 103C, 503C & 504D, Technology Building & No. 3A &
4A, Technology Building Annex, Zhaoshang Sub-District,
Nanshan District,
Shenzhen, China**

***in vitro* diagnostic medical devices
self-testing**

SARS-CoV-2 Antigen Kit (Colloidal Gold)

REF: CG123001, CG123003, CG123005, CG123007, CG123025

**in terms of design documentation, comply with requirements
of Annex III (Section 6) to Directive 98/79/EC (as amended)
implemented into Polish law,
as evidenced by the audit conducted by the PCBC**

Validity of the Certificate: from 10.11.2021 to 27.05.2024

The date of issue of the Certificate: 10.11.2021

The date of the first issue of the Certificate: 10.11.2021



Issued under the Contract No. **MD-40/2021**
Application No: **397/2020**
Certificate bears the qualified signature.
Warsaw, 10/11/2021
Module A1

**Anna
Małgorzata
Wyroba**

**Vice-President
Mgr. Anna Wyroba**

**Elektronicznie
podpisany przez Anna
Małgorzata Wyroba
Data: 2021.11.10
16:14:07 +01'00'**

Other 3rd party certificates and recommendations

EU common list

<https://covid-19-diagnostics.jrc.ec.europa.eu/devices/detail/1197>



EUROPEAN COMMISSION
DIRECTORATE-GENERAL FOR HEALTH AND FOOD SAFETY
Public health, country knowledge, crisis management
Health Security

EU health preparedness:

A common list of COVID-19 rapid antigen tests;
A common standardised set of data to be included in COVID-19 test result certificates; and
A common list of COVID-19 laboratory based antigenic assays

Agreed by the Health Security Committee

Manufacturer	RAT commercial name	Device ID # ¹¹	Clinical performance As reported by independent validation studies	Clinical performance Data by manufacturer	EU Member States using in practice	Other countries using in practice	Completed validation studies	SARS-CoV-2 Target protein	Included in EU common list since:
Goldsite Diagnostic Inc.	SARS-CoV-2 Antigen Kit (Colloidal Gold)	1197	<i>Retrospective in vitro study</i>						
			DE: Positive evaluation by Paul-Ehrlich-Institut: Sensitivity of 100% at Ct ≤ 25; Manufacturer specificity: 100%	93.04% sensitivity; 100% specificity Nasal swab	BE, BG, CY, FR, RO, SI, ES	UK	FR, DE ¹² , ES UK	Unknown	14 July 2021

COVID-19 In Vitro Diagnostic Medical Device - detail

SARS-CoV-2 Antigen Kit (Colloidal Gold)

Manufactured by Goldsite Diagnostics Inc, China - en.goldsite.com.cn/ 

Device identification number	1197
CE Marking	 Yes
HSC common list	 Yes
HSC mutual recognition	 Yes
Format	Near POC / POC
Target	Antigen
Specimen	Nasal swab, Other
Commercial Status	Commercialised
Last Update	2021-07-13 06:05:32 CET
Comments	

Hide HSC list status history 

[View current HSC list](#)

2021-07-13 Added to mutual recognition list

Germany Paul-Ehrlich-Institut Evaluation

<https://www.pei.de/EN/newsroom/dossier/coronavirus/test-systems.html>

Bundesinstitut für Impfstoffe und biomedizinische Arzneimittel
Federal Institute for Vaccines and Biomedicines

Paul-Ehrlich-Institut 

25.06.2021

Comparative evaluation of the sensitivities of SARS-CoV-2 antigen rapid tests

Goldsite COVID-19 SARS-CoV-2 Antigen Kit (Colloidal Gold)	Goldsite Diagnostics Inc.
--	---------------------------

France ANSM self-test special approval

<https://covid-19.sante.gouv.fr/tests>



MINISTÈRE
DES SOLIDARITÉS
ET DE LA SANTÉ

PLATEFORME COVID-19

Se connecter



VISUALISATION TESTS COVID-19

Nombre de tests

1/518

Les anciens kits de criblage 1^{er} et 2^e ligne qui ne
diffèrent pas les 3 nouvelles mutations conformément
au MINSANTE N°2021-72 du 26 mai, seront retirés de
la liste à compter du 27 juin 2021

Cliquez pour
accéder à la liste
commune
européenne TAB

Signalement

Problème

Contextes juridiques

Cliquez pour sélectionner et télécharger les fichiers des
contextes juridiques

Statut

☐ CE ☐ CMR ☐ ME ☐ HES

Type de test

Sous-type de test

Code

Type prélevement

Rechercher

Rechercher

Tableau de bord des tests

Cliquez pour déplacer et modifier les profils du
tableau de bord

1 test affiché

Options

NOM	FABRICANT	DISTRIBUTEUR	CE	UE	CMR	SOUS-TYPE DE TEST	CODE	TYPE DE PRÉLEVEMENT
SWAB-CAR-0-Kit antigénique-DV Cathedra	Goldsite Diagnostics Inc		☐	☐	☐	Autocrit	fr	Saliva

Belgian recommendation list

https://www.afmps.be/fr/humain/produits_de_sante/dispositifs_medicaux/covid_19/tests



Manufacturer	Test	Sensitivity	Specificity	Specimen type	Test ID in JRC database	Test on common HSC list?
Goldsite Diagnostics	SARS-CoV-2 Antigen Kit (Colloidal Gold)	93.0	100.0	Nasal swab	1197	Yes

UK DHSC Phase 3A validation

<https://www.gov.uk/government/publications/assessment-and-procurement-of-coronavirus-covid-19-tests/outcome-of-the-evaluation-of-rapid-diagnostic-assays-for-specific-sars-cov-2-antigens-lateral-flow-devices>



Department
of Health &
Social Care

Guidance

Outcome of the evaluation of rapid diagnostic assays for specific SARS-CoV-2 antigens (lateral flow devices)

Updated 11 October 2021

Table 1: summary of lateral flow devices that have passed phase 3a validation

Lateral flow device	Status	Date evaluation completed
Goldsite SARS-CoV-2 Antigen Kit (Colloidal Gold)	Pass	8 July 2021



GOLDSITE DIAGNOSTICS INC.
No. 103C, 503C & 504D, Technology Building &
No. 3A & 4A, Technology Building Annex,
Zhaoshang Sub-District, Nanshan District,
Shenzhen, China, 518067
Email: export@goldsite.com.cn
Website: www.goldsite.com.cn

Malaysia recommendation list

<https://portal.mda.gov.my/announcement/596-list-of-recommended-for-use-of-covid-19-ivd-test-kit.html>

List of Recommended for Use of Covid-19 IVD Test kit (Professional use only)

SPECIAL ACCESS LIST OF COVID-19 TEST KIT (FOR PROFESSIONAL USE ONLY)

The list of Covid-19 Antibody, Antigen and Nucleic Acid Amplification Test (NAAT) that is recommended for use based on the decision on the consensus of the Covid-19 Test Kit Expert committee is as follows:

All test submissions are scored according to:-

- the manufacturer reported clinical and analytical performance
- the evaluation test results from testing facilities are according to the committee evaluation criteria set by Clinical expert
- Supporting Documents for COVID-19 IVD Test Kits **Special Access**

Below is the list of all tests that have been selected to date and the status is Recommended for Use (please note: list is updated on a routine basis)

Notes:

*Sample type and lot number is based on testing facility evaluation report.

1%	Vega Gemilang Enterprise	Goldsite SARS-CoV-2 Antigen kit	Goldsite Diagnostic Inc. China	CG123025	20210415	RTK Antigen	Nasal Swab
----	--------------------------	--	---------------------------------------	----------	----------	-------------	------------

Thai FDA self-test approval

https://www.fda.moph.go.th/sites/Medical/SitePages/test_kit_covid19.aspx



กองควบคุมเครื่องมือแพทย์

Medical Device Control Division

สำนักงานคณะกรรมการอาหารและยา

 รายชื่อชุดตรวจสำหรับ COVID-19 ประเภท Rapid Test Antigen หรือ Antigen Test Kits แบบตรวจแอนติเจนด้วยตนเอง (COVID-19 Antigen Test Self-Test Kits) ที่ได้รับการอนุญาตให้ผลิตนำเข้า จากสำนักงานคณะกรรมการอาหารและยา ข้อมูล ณ วันที่ 22 พฤศจิกายน 2564						
ลำดับที่	ชื่อผลิตภัณฑ์	ชื่อผู้นำเข้า	ชื่อผู้ผลิต	วันที่ได้รับ อนุญาต (วัน/เดือน/ปี)	เลขที่ใบรับขอ ประเมิน เภส บาล	ลิงก์วิดีโอ ขั้นตอนการใช้ ชุดตรวจ
138	SARS-CoV-2 Antigen Kit (Colloidal Gold) รหัสสินค้า CG123001, CG123005 ขนาดบรรจุ 1, 5 ชุดทดสอบแต่ละชุด (Vial) (swab) 	จำเริญส่วนจำกัด โอ โย แอตาบิล โทร. 0 3428 7410	GOLDSITE DIAGNOSTICS INC. China	18/11/2564	T 6400536	 



DECLARATION OF CONFORMITY

Manufacturer	Goldsite Diagnostics Inc.	
Address	No.103C, 503C & 504D, Technology Building & No. 3A & 4A, Technology Building Annex, Zhaoshang Sub-District, Nanshan District, Shenzhen, China, 518067	
European Representative	CMC MEDICAL DEVICES & DRUGS, S.L. C/ Horacio Lengo No 18, CP 29006, Málaga-Spain	
Product Information	SARS-CoV-2 Antigen Kit (Colloidal Gold) (for self-testing)	
GMDN terms	SARS-CoV-2 antigen IVD, kit, immunochromatographic test (ICT), self-testing. GMDN Code: 65454	
Conformity Assessment Route: Annex III, point 6	<i>We herewith declare that the above-mentioned products meet the provisions of the following EC Council Directives and Standards. All supporting documentations are retained under the premises of the manufacturer.</i>	
General Applicable Directives:	<i>In Vitro Diagnostic Medical Devices DIRECTIVE 98/79/EC OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL Of 27 October 1998 Classification: Device for self-testing Directive 98/79/EC, Article 9, Annex III</i>	
Standards Applied	EN ISO 13485:2016 BS EN 13612:2002 ISO 14971:2019 EN ISO 23640:2015	ISO 15223-1:2016 EN ISO 18113-1: 2011 EN ISO 18113-4: 2011 EN 62366-1: 2015



Place, date of issue: Shenzhen, P.R. China, November 10, 2021

Signature of General Director

SIGNATURE



STAMP

BOX Label

1T


REF CG123001




BOX Label

5T

REF CG123005





SARS-CoV-2 Antigen Kit (Colloidal Gold) (for self-testing)

Intended Use

This SARS-CoV-2 Antigen Kit (Colloidal Gold) is a single use, in vitro visually read rapid immunochromatographic assay for the qualitative detection of SARS-CoV-2 nucleocapsid (N) antigen in human nasal swab specimens. Results are for the identification of SARS-CoV-2 nucleocapsid protein antigen. Antigen is generally detectable in nasal swabs during the acute phase of infection. Positive results indicate the presence of viral antigens, but correlation with patient history and other diagnostic information is necessary to determine the infection status. Positive results do not rule out bacterial infection or co-infection with other viruses. The agent detected may not be the definite cause of disease.

Negative results should be treated as presumptive and confirmation with a molecular assay, if necessary, for patient management may be performed. Negative results do not rule out COVID-19 and should not be used as the sole basis for treatment or patient management decisions, including infection control decisions. Negative results should be considered in the context of patient's recent exposures, history and the presence of clinical signs and symptoms consistent with COVID-19. Persons who test negative and continue to experience COVID-19 like symptoms of fever, cough and/or shortness of breath may still have SARS-CoV-2 infection and should seek follow up care with their physician or healthcare provider.

The SARS-CoV-2 Antigen Kit (Colloidal Gold) is intended to be used manually by lay users (self-testing) in a private setting to aid in the rapid diagnosis of SARS-CoV-2 infections. Children (at least 2 years old) younger than 15 years old should be tested by an adult. Rapid diagnosis of SARS-CoV-2 infection will help healthcare professionals to treat patients and control the disease more efficiently and effectively.

Summary and Explanation

The novel coronaviruses belong to the β genus. COVID-19 is an acute respiratory infectious disease. People are generally susceptible. Currently, the patients infected by the novel coronavirus are the main source of infection; asymptomatic infected people can also be an infectious source. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. The main manifestations include fever, fatigue and dry cough. Nasal congestion, runny nose, sore throat, myalgia and diarrhea are found in a few cases.

Principle of the Procedure

The kit employs lateral flow immunoassay technology. Using this test allows for the rapid detection of nucleocapsid protein from SARS-CoV-2. To begin the test, a self-collected nasal swab sample is inserted into the Extraction Tube. The pre-filled solution in the tube interacts with the sample and facilitates exposure of the appropriate viral antigens to the antibodies used in the test. The extracted sample is then added to the sample well of the Test Cartridge.

If the extracted sample contains SARS-CoV-2 antigens, a pink-to-red line will appear on the test line region (T) and a red line will appear on the procedural control line region (C). Absence of the test line (T) indicates SARS-CoV-2 antigens are not detected, suggesting SARS-CoV-2 is not present, or is present at very low levels.

Materials Required But Not Provided

1. Timer
2. Hand soap and water or hand sanitizer
3. Household waste container

Materials Provided

Catalog number	Test Cartridge	Sampling Swab	Extraction Tube	Instructions for Use
CG123001	1	1	1	1
CG123003	3	3	3	1
CG123005	5	5	5	1
CG123007	7	7	7	1
CG123025	25	25	25	1

Warnings and Precautions

- For in vitro diagnostic use.
- Read the instructions fully before starting the test. To obtain accurate results, the instructions must be followed.
- If the package has been damaged, the label cannot be seen clearly or if the cartridge has expired, do not use the cartridge.
- Do not open the foil pouch of the Test Cartridge until ready for use.
- Do not eat the desiccant.
- Do not reuse the cartridges, tubes or swabs.
- Do not interchange or mix components from different kit lots.
- Do not touch swab tip when handling the swab.
- Inadequate or inappropriate specimen collection, may yield false negative test results.
- Dispose of kit components and samples in household trash.
- Use of gloves during sampling and testing is recommended.

Storage and Stability

Store the kits at 2 – 30° C in a dry place and avoid direct sunlight. Do not freeze any of the kit components. Keep the kits out of reach of children. The unopened cartridges are stable until the expiry date printed on the labels. Once the foil pouch of the test cartridge is opened, the test should be started within 30 minutes.

Test Procedure

Note: Before starting the test, wash your hands with soap and water or use hand sanitizer. Make sure they are dry before starting.

1. Allow the kit to equilibrate to room temperature (15 – 30°C) before testing. Open the kit and identify kit components and instructions.

2. Collecting the nasal swab sample. Open swab package at stick end. Take swab out. While gently rotating, insert the swab into one nostril. The swab tip should be inserted up to 2 cm (1/2 to 3/4 of an inch) from the edge of the nostril. Rotate the swab 5 times against the nasal wall. Remove and repeat the sampling process using the same swab for the other nostril.

Note: Do not touch swab tip when handling the swab.
Note: Children (at least 2 years old) younger than 15 years old, and people who are unable to perform the test themselves including the elderly and the sick should be tested by another adult. To sample a child, insert the swab into one of their nostrils until you feel some resistance. Rotate the swab 5 times against the nasal wall. Remove the swab and insert the same swab into the other nostril, repeat the sampling process. Do not continue the test if the child feels any pain.



Note: Inadequate or inappropriate sample collection may yield false negative results.

3. Remove the Dropper Lid from the purple pre-filled Extraction Tube. Place and soak the swab into the tube.

Remove the Dropper Lid

Place and soak the swab into the tube

4. Rotate the swab 5 – 6 times while squeezing the sides of the tube. Insert the tube into the hole indicated on the kit box. Make sure the tube is standing upright and reaches the bottom. Leave the swab in the Extraction Tube for 1 minute.

Rotate the swab 5 – 6 times while squeezing

5 – 6 times Squeeze

Insert the tube to the hole on the kit box. Leave the swab for 1 minute.




5. Stir to mix the contents well. Remove the swab while squeezing the sides of the tube. Immediately discard the swab into the garbage.

Stir to mix the contents well

Stir

Remove the swab while squeezing the sides of the tube



6. Attach and tighten the Dropper Lid to the top of the Extraction Tube. Mix the contents thoroughly (by inverting the tube several times).

Attach and tighten the Dropper Lid

Invert gently several times



7. Open the foil pouch of the Test Cartridge. Place the cartridge on a clean, flat and dry surface. Label the cartridge with patient ID number. Remove the cap on top of the Dropper Lid, invert the extraction tube, and then add two drops (around 70 uL) of the well-mixed sample into the sample well of the cartridge.

Remove the white cap

Add two drops into sample well




8. Leave the sample-loaded cartridge at room temperature for 12 minutes. Do not handle or move the cartridge during this time.

9. After the 12-minute incubation, read the results. Do not interpret the results after 15 minutes (from addition of the sample).

Interpretation of Results

Positive Result:

Two lines will appear.
One pink-to-red line on the test line region (T) and one red line on the control line region (C).
Note: Look very closely! The test line (T) can be very faint. Any pink/red visible T line should be considered positive.



Negative Result:

A single line appear on the control line region (C).



Invalid Result:

No red line appears on the control line region (C).
Note: Insufficient sample volume, incorrect operating procedure are the most likely reasons for control line failure. If the test result is invalid, a new swab should be collected, and the test should be performed again with a new kit. If the problem persists, discontinue using the kit immediately and contact your local distributor.



Limitations

- The test is only for qualitative detection of the SARS-CoV-2 antigen in nasal swab samples.
- This test detects both viable (live) and non-viable SARS-CoV-2. Test performance depends on the amount of virus (antigen) in the sample and may or may not correlate with viral culture results performed on the same sample.
- A negative test result may occur if the level of antigen in a sample is below the detection limit of the test.
- Failure to follow the Test Procedure, Interpretation of Results may adversely affect test performance and/or invalidate the test results.
- False negative results may occur if a sample is improperly collected or handled.
- Positive test results do not rule out co-infections with other pathogens.
- A negative result does not exclude SARS-CoV-2 infection. Negative results should be treated as presumptive and may need to be confirmed by a molecular assay.
- False negative results are more likely after eight days of symptoms.

Performance

Clinical Performance

The clinical performance of SARS-CoV-2 Antigen Kit (Colloidal Gold) (for self-testing) was determined by testing samples from 267 persons suspected of COVID-19. The samples were collected within 7 days post-onset of symptoms or suspected exposure. For each individual, two swabs were collected. The first one was self-sampled nasal swab and was self-tested directly with SARS-CoV-2 Antigen Kit (Colloidal Gold) (for self-testing) at the site. The second was a healthcare-collected nasopharyngeal swab which was placed in virus transport medium, shipped to laboratory, and determined to be positive or negative using an NMFA (National Medical Products Administration, China) approved RT-PCR method, i.e., the comparator method.

		RT-PCR Test Results		
		Positive	Negative	Total
SARS-CoV-2 Antigen Kit (Colloidal Gold)	Positive	110	1	111
	Negative	9	147	156
	Total	119	148	267

Clinical sensitivity: 92.44% (95% CI: 86.13 – 96.48%)
Clinical specificity: 99.32% (95% CI: 96.29 – 99.98%)
Overall percent agreement: 96.25% (95% CI: 93.22 – 98.19%)

Limit of Detection (LoD)

SARS-CoV-2 Antigen Kit (Colloidal Gold) was confirmed to detect 2.5 ng/mL of SARS-CoV-2 nucleocapsid protein antigen.

Cross Reactivity and Microbial Interference

The cross-reactivity with the following microorganisms was examined. Samples that tested positive for the following microorganisms were negative when tested with the SARS-CoV-2 Antigen Kit (Colloidal Gold). The microbial interference study evaluated whether microorganisms possibly contained in clinical samples interfere with the detection capability of the kit which may lead to false negative results. Each microorganism was tested in triplicate in the presence of a fabricated SARS-CoV-2 positive sample (concentration: $3 \times \text{LOD}$). No cross-reactivity or interference with the microorganisms listed in the table below was found.

No.	Microorganism	Final Test Concentration
1	HCcV-CC43	$2 \times 10^6 \text{ TCID}_{50}/\text{mL}$
2	HCcV-229E	$2 \times 10^6 \text{ TCID}_{50}/\text{mL}$
3	HCcV-NL63	$2 \times 10^6 \text{ TCID}_{50}/\text{mL}$
4	RSV	$2 \times 10^6 \text{ TCID}_{50}/\text{mL}$
5	Rotavirus	$2 \times 10^6 \text{ TCID}_{50}/\text{mL}$
6	MERS	$1 \times 10^6 \text{ TCID}_{50}/\text{mL}$
7	Adenovirus	$2 \times 10^6 \text{ TCID}_{50}/\text{mL}$
8	Norovirus	$2 \times 10^6 \text{ TCID}_{50}/\text{mL}$
9	Mycoplasma pneumonia	$1.5 \times 10^8 \text{ cfu}/\text{mL}$
10	Influenza A Virus (H1N1)	$2 \times 10^7 \text{ TCID}_{50}/\text{mL}$
11	Influenza B virus (Yamagata)	$2 \times 10^7 \text{ TCID}_{50}/\text{mL}$

Interference

The following interfering substances have no impact on SARS-CoV-2 Antigen Kit (Colloidal Gold).

No.	Interfering substance	Final Test Concentration
1	Phenylephrine	15% v/v
2	Oxymetazoline	15% v/v
3	Chloride de Sodium	5 mg/mL
4	Bedomethasone	5 ng/mL
5	Dexamethasone	0.5 µg/mL
6	Fluticadide	0.5 µg/mL
7	Triamcinolone acetonide	1 ng/mL
8	Eudesonide	2.5 ng/mL
9	Mometasone	1 ng/mL
10	Fluticasone	2 ng/mL

Hook Effect

There is no hook effect at 600 µg/mL of SARS-CoV-2 nucleocapsid protein antigen.

Frequently Asked Questions

1. Will this test hurt?

No, the nasal swab is not sharp and it should not hurt. Sometimes the swab can feel slightly uncomfortable. If you feel pain, please stop the test and seek advice from a healthcare provider.

2. What are the known and potential risks and benefits of this test?

Potential risks include:

- Possible discomfort during sample collection.
- Possible incorrect test results.

Potential benefits include:

- The results, along with other information, can help your healthcare provider make informed recommendations about your care.
- The results of this test may help limit the spread of COVID-19 to your family and others in your community.

3. What does it mean if I have a positive test result?

If you have a positive test result, it is very likely that you have COVID-19 because proteins from the virus that causes COVID-19 were found in your sample. Therefore, it is also likely that you may be placed in isolation to avoid spreading the virus to others. There is a very small chance that this test can give a positive result that is wrong (a false positive result). If you test positive with the SARS-CoV-2 Antigen Kit (Colloidal Gold) you should self-isolate and seek follow-up care with your healthcare provider as additional testing may be necessary. Your healthcare provider will work with you to determine how best to care for you based on your test result(s) along with your medical history, and your symptoms.

4. What does it mean if I have a negative test result?

A negative test result means that proteins from the virus that causes COVID-19 were not found in your sample. It is possible for this test to give a negative result that is incorrect (false negative) in some people with COVID-19. This means that you could possibly still have COVID-19 even though the test is negative. The amount of antigen in a sample may decrease the longer you have symptoms of infection. Specimens collected after you have had symptoms for more than five days may be more likely to be negative compared to a molecular assay.

If you test negative and continue to experience COVID-19 like symptoms of fever, cough and/or shortness of breath you should seek follow up care with your healthcare provider. For example, your healthcare provider may suggest you need another test to determine if you have contracted the virus causing COVID-19. It is important that you work with your healthcare provider to help you understand the next steps you should take.

5. What is the difference between an antigen and molecular test?

An antigen test, such as the SARS-CoV-2 Antigen Kit (Colloidal Gold), detects proteins from the virus. Molecular tests detect genetic material from the virus. Antigen tests are very specific for the virus, but not as sensitive as molecular tests. This means that a positive result is highly accurate, but a negative result does not rule out infection. If your test result is negative, you should discuss with your healthcare provider on whether an additional test is necessary and if you should continue isolating at home.

Symbols

	Consult instructions for use		Tests per kit		Authorized representative in the European Community
	In vitro diagnostic medical device		Use-by date		Do not reuse
	Temperature limit 2 - 30°C		Lot number		Catalog number
	Avoid sunshine		Manufacture		Date of manufacture

Goldsite Diagnostics Inc.

Address of manufacturer
No. 103C, 503C & 504D, Technology Building & No. 3A & 4A, Technology Building Annex, Zhaoshang Sub-District, Nanshan District, Shenzhen, China, 518067
Manufacturing Site
No. 103C Technology Building & No. 3A & 4A, Technology Building Annex, Zhaoshang Sub-District, Nanshan District, Shenzhen, China, 518067
Tel: 86 755 26890807
Fax: 86 755 26890799

CEC MEDICAL DEVICES & DRUGS, S.L.
C/ Horacio Lango No 18, CP 29006, Málaga Spain

CE 1434

IFU4451601
Version: V2
16/11/2021